

# BiAA Upgrade Instruction

# BiAA Introduction

## BiAA – Built-in AEC Assistance

With BiAA, AEC exposures can also be performed for non-Bucky imaging, eliminating the need for manual exposure settings and thereby reducing the risk of retakes and extra patient dose.

For Bucky imaging the built-in AEC chamber is still used.

BiAA is available with the following detectors:

- CXDI-420C Wireless
- CXDI-720C Wireless
- CXDI-820C Wireless



# BiAA – *System Requirements*

Basic Requirements for upgrade:

- Arcoma system with Canon Elite detector & Multibox MB-02.
- Canon NE 3.10 SW or later > *Shall be upgraded to version according to Arcoma Software Release note*
- Latest Arcoma SW (see Arcoma Partner portal) and corresponding / compatible Softwares Canon and CPI.



# BiAA Upgrade: Preparation

## 1) Download from ARCOMA Partner portal

- a) BiAA Upgrade Instruction (*this*)
- b) BiAA package (Software, xml-file, text-file)
- c) System Software for Latest Version Upgrade

## 2) Collect information from hospital + Update SW [→ Info](#)

### > Hospital Visit or by Remote Access

- a) Upgrade system with Canon *BuiltInAECAssistance* SW and Generate a Canon License Request File with BiAA option selected.
- b) Create a copy of the CANONKIT.txt-file.
- c) Check CPI generator serial number and System ID.

## 3) Place an order to Arcoma

Article number: 0073-815-011 Upgrade BiAA

We need the following information from you:

- a) Serial number of the system
- b) Serial number and System ID of the generator
- c) CANONKIT.txt-file
- d) Canon License Request File (xx.lrf-file)

## 4) Arcoma will delivery to you:

- a) CPI License key
- b) CANONKIT.lic-file
- c) Canon License key



# BiAA Upgrade: Perform

**4) Upgrade the System SW to the latest released version.**

**5) Perform the BiAA Upgrade** → [Info](#)

- a) CANONKIT.lic-file update (CXDI\_NE\_Overwrap folder)
- b) Settings in GenConfig (GenConfig)
- c) Generator license key update (Genware)
- d) Receptor Settings (Genware)
- e) Register the Canon license key (Canon Service Tool)
- f) Add and update files in CCS-S-folder:
  - a) DRTSETTING.ini.file
  - b) CanonWarnings and safety check xml files
  - c) BiAA Error Message update

**6) Adjust BiAA cut-off dose** → [Info](#)

**7) BiAA settings in Anatomical Protocols** → [Info](#)





# Collect Information for BiAA License Update

# Required Information for BiAA License Update

**The listed information is required for BiAA License Update. See following pages on how to collect the required information.**

- a) Serial number of the system
- b) Serial number and System ID of the generator
- c) CANONKIT.txt-file
- d) Canon License Request File (xx.lrf-file)

Information can be collected via Hospital visit or by Remote Access to the system.



# System serial number and Generator Serial Number


## SYSTEM SERIAL NUMBER

> See labels on the System Cabinet.

## GENERATOR SERIAL NUMBER

**WARNING!** *Read Warnings in System Manual before opening up the System Cabinet.*

- Open the System cabinet covers,  
> *Read the Manual!*
- Generator Serial Number is found on the right side of the generator.



Generator Serial Number





# System ID



## SYSTEM ID

1. Open GenWare
2. Select LicKey button
3. System ID = numbers in the green fields

Generator License Key Setup

System ID / License Key

System ID 0000 - 0000 - 1421 - 383a

License Key | - - -

Apply

Close

0 1 2 3 4 5 6 7 8 9

Q W E R T Y U I O P BACKSPACE


A S D F G H J K L CLEAR

Z X C V B N M Lower Case Letters

# CANONKIT.TXT-file

## CANONKIT.TXT-file

1. Go to CXDI\_NE\_Overwrap folder
2. Copy: CANONKIT.TXT-file



# Canon License Request File

- 1) Install the BuiltInAECAssistance SW
  - a) Start the Setup.exe file
  - b) Select Generator Console Mode as Built-in AEC Configuration
  - c) Follow the steps in the installation tool until finished.
- 2) Open Canon Service Tool.
- 3) Canon Service Tool: **Add-in Settings** / Built-in AEC Assistance  
> Select *Enable settings of ROI and exposure*.
- 4) Canon Service Tool: **License**. Select *Create License Request Files* and select *Option Software Built-in AEC Assistance* > Confirm.
- 5) A lrf-file is created.

The screenshot illustrates the steps to create a license request file in the Canon Service Tool. It shows three windows: 'Menu selection', 'Add-in software setting', and 'CREATE LICENSE REQUEST FILES'.

**Menu selection:** The 'Add-in Setting' and 'License' icons are highlighted with a yellow box. A red arrow points from this box to the 'Add-in software setting' window.

**Add-in software setting:** The 'Built-in AEC Assistance' tab is selected. The 'Module setting' section shows 'Enable settings of ROI and exposure' checked. A red arrow points from this box to the 'CREATE LICENSE REQUEST FILES' window.

**CREATE LICENSE REQUEST FILES:** The form is filled out with the following information:

- Institution Name: Arcoma AB
- Installed Room: Demo room Precision i5
- Institution Department:
- Street Address:
- City:
- State/Province/Prefecture:
- Country or Region: Sweden
- Postal Code (ZIP): SE-352 46 Vaxjo
- Dealer (ID/NAME): Arcoma AB
- License Request Comment:

Below the form, the 'Select the purchased optional software' section shows the following options:

System Software	Option Software
<input checked="" type="checkbox"/> CXDI Control Software NE	<input checked="" type="checkbox"/> Scatter Correction for CXDI Series
	<input checked="" type="checkbox"/> Advanced Edge Enhancement
	<input checked="" type="checkbox"/> Free Rotation for CXDI Series
	<input checked="" type="checkbox"/> Built-in AEC Assistance

A red arrow points from the 'Option Software' section to the 'Confirm' button.

# Perform the BiAA Upgrade

# BiAA Upgrade – Steps

**Upgrade the System SW to the latest released version.**

## **5) Perform the BiAA Upgrade**

- a) CANONKIT.lic-file update (CXDI\_NE\_Overwrap folder)
- b) Settings in GenConfig (GenConfig)
- c) Generator License key update (Genware)
- d) Receptor Settings (Genware)
- e) Register the Canon license key (Canon Service Tool)
- f) Add and update files in CCS-S-folder:
  - a) DRTSETTING.ini.file
  - b) CanonWarnings and safety check xml files
  - c) BiAA Error Message update






# CANONKIT.LIC-file update

## CANONKIT.lic-file UPDATE



1. Go to CXDI\_NE\_Overwrap
2. Exchange the current CANONKIT.LIC file with the new provided file.



# GenConfig Settings

## GENCONFIG SETTINGS UPDATE

1. Open GenConfig (CXDI\_NE\_Overwrap folder)
2. Select tab: Licensed Features
3. Select 5-Field BiAA



# CPI License Key



## License Key Update

1. Open GenWare
2. Select LicKey button
3. Enter the provided License Key

Generator License Key Setup

System ID / License Key

System ID: 0000 - 0000 - 1421 - 383a

License Key: | - - -

Apply

Close

0 1 2 3 4 5 6 7 8 9

Q W E R T Y U I O P BACKSPACE


A S D F G H J K L CLEAR

Z X C V B N M Lower Case Letters

# CPI License Key

## License Key Update

1. Open GenWare
2. Select LicKey button
3. Enter the provided License Key
4. Check that option 3 Flat Panel AEC is shown




# GENWARE, Receptor

AEC Channel = 5 shall be selected when BiAA used.

## STANDARD SETUP

- R1: Table
- R2: Table top
- R3: Wall stand
- R4: Wall stand detector, out of Bucky
- R5: Free detector

R2, R4 and R5 recommended for BiAA use.






# CANON LICENSE UPDATE

## Canon BiAA License Update


1. Open Canon Service Tool
2. Canon Service Tool: License.  
Select *Register license keys* and enter the provided license key.
3. Check that option Built-in AEC Assistance is shown as Activated.



# Add and update files in CCS-S folder

**Add / Update the following files in the CCS-S folder:**

- DRTSETTING.ini-file
- UsePixelValueAlertWithBiAA.xml
- UseWifiAlertWithBiAA.xml
- Exchange the following file in the CCS-S/en-US folder:  
*Canon.Medical.DR.Seychelle.Facade.Common.Properties.ErrorMasgResources.en-US.resources*



# BiAA Settings

# Canon APR Editor: BiAA settings

Anatomical Protocols are defined in the Canon APR Editor.

## Settings required to use BiAA in a protocol:

- Technique: AEC
- ms = Used as Backup value
- Detector AEC Assist: YES

## Adjustable protocol settings related to BiAA:

- Active Detector AEC
- DRT Density
- Detector Rotation

*See next page for more information.*

Used as  
Backup value

Technique: AEC

Parameters		Configuration			
	NAME	Very Small	Small	Medium	Large
▶	Rad kV	40	40	95	40
	Rad mA	320.0	250.0	320.0	320.0
	ms	16.0	16.0	16.0	16.0
	mAs	5.1	4.0	5.1	5.1
	Technique	AEC	AEC	AEC	AEC
	Film	Film Screen 1	Film Screen 1	Film Screen 1	Film Screen 1
	Focus	LARGE	LARGE	LARGE	LARGE
	Left Field	NO	NO	NO	NO
	Center Field	NO	NO	NO	NO
	Right Field	NO	NO	NO	NO
	Receptor	4	4	4	4
	Density	0	0	0	0

New parameters  
related to BiAA

	Detector AEC Assist	YES	YES	YES	YES
	Active Detector AEC	A	A	A	A
	DRT Density	0	0	0	0
	Detector Rotation	0 degree	0 degree	0 degree	0 degree

# BiAA – Protocol Setup

Parameters<sup>\*)</sup> related to BiAA in Anatomical Protocol:

	Detector AEC Assist	YES
	Active Detector AEC	C
	DRT Density	0
	Detector Rotation	0 degree

**Shall be YES to activate BiAA**

*> Can be changed by the user*

**Select active chambers**




*> Can be changed by the user*

**Adjust the cut-off dose**

*> Impact on image noise level*

**BiAA mode:**

**Manual/0 degree or Auto Rotation link**



<sup>\*)</sup> Note that there are also other parameters related to BiAA in the standard configuration, recommendation to hide these parameters with the support of the Configuration tab settings.



# Define AEC cut-off dose

## Check the BiAA cut-off dose:

- 1) Select a BiAA protocol and activate the center chamber ( C ).
- 2) Position the X-ray tube above the detector (SID 115 cm) and adjust the collimator light to cover the detector area.
- 3) Make an exposure with the following parameters  
- RQA5: 70 kV, 21 mmAl added filtration
- 4) Note the EI value received (shown in Canon, exposed image).

*Example of EI-values:*

EI : 200 = 2  $\mu$ Gy cut-off dose

EI : 180 = 1,8  $\mu$ Gy cut-off dose

$EI = 100 \times \text{Base} / \text{Sfpd}$ ,  $\text{Sfpd} = 678 \text{ LSB}/\mu\text{Gy}$  at RQA5

## Adjust the BiAA cut-off dose:

- 1) Open the drtsetting.ini-file in the Canon CCS folder.
- 2) Adjust the Base value until correct cut-off dose received.

Base = 1350 > EI = 200 > cut-off dose  $\sim 2\mu\text{Gy}$

Base = 2115 > EI  $\sim 312$  > cut-off dose  $\sim 3\mu\text{Gy}$



# DRTSETTING.ini-file

The drtsetting.ini-file shall be saved to Canon CCS folder.



Parameters of interest in the drtsetting:

- **Base** > Defines the AEC cut-off dose
- **Early Termination Parameters** > Used to terminat the exposure if it is not reaching the selected AEC chamber.
- **Density** > Possibility to adjust the AEC cut-off dose in each Protocol.
- **Sensitivity** > Adjustment of the cut-off dose

Range  
Settings

```
[COMMON]
SectionNum=2
BaseRange=0
SensRange=0
DensRange=4
EarlyTermination=1
ETThreshold=4
ETTime=20
```

Not used

BiAA with the B1 detector  
(not implemented)

```
[DRT0]
Base=1
BaseMax=100
SensMax=100
DensMax=100
```

```
[DRT1]
```

```
Base=1350
```

Max value  
Settings

```
BaseMax=200
SensMax=200
DensMax=200
```

Time  
settings

```
DirectionAuto=1
SystemDelay=1000
SystemDelay_Wireless=1200
ETStartTime=1100
ETStartTime_Wireless=1300
```

**Do not change!**

Used to define selection  
of AEC chambers in the  
Protocol editor.

```
[LUT1]
```

```
C=01,02,03,04,05,10,11,12,13,14,15,16,17,18,19,20,21,22,23,24,25,
F=0004,0001,0010,0100,0040,0005,0140,0014,0011,0110,0050,0104,004
```

# THANK YOU

Arcoma AB | Annavägen 1 | 352 46 Växjö | Sweden



## PHONE

+46 470 70 69 00



## EMAIL

[office@arcoma.se](mailto:office@arcoma.se)



## WEBSITE

[www.arcoma.se](http://www.arcoma.se)

# DISCLAIMER

The information presented in this training material is intended solely for educational and informational purposes. The content provided is based on the knowledge and expertise available up to the date of this presentation. While every effort has been made to ensure the accuracy and completeness of the information contained herein, no representation or warranty, expressed or implied, is made regarding the accuracy, reliability, suitability, or completeness of the content.

The material provided does not constitute professional advice, whether legal, financial, medical, or otherwise. It is recommended that you consult with appropriate professionals or experts in the relevant field before making decisions or taking actions based on the information provided in this training material.

Any reliance you place on the material presented in this training material is strictly at your own risk. The creators, authors, and presenters of this material disclaim any liability for any loss, damage, or injury arising from the use or misuse of the information provided. Furthermore, the views and opinions expressed in this training material are those of the authors and presenters and do not necessarily reflect the official policy or position of any organization, company, or institution.

By using this training material, you acknowledge and agree to the terms of this disclaimer. If you do not agree with these terms, please refrain from using this material.

**Copyright:** All rights reserved. No part of this training material may be reproduced, distributed, or transmitted in any form or by any means, including photocopying, recording, or other electronic or mechanical methods, without the prior written permission of the creators. Please consult the appropriate authorities or legal advisors if you have any questions regarding the use or distribution of this training material.