


WELCOME TO TODAY'S WEBINAR FROM



WORLD CLASS DIGITAL X-RAY SYSTEMS


The presentation will soon begin. Please prepare by disabling your camera and microphone:



VIDEO MUTE
OFF

THE PRACTICAL STUFF:

- VIDEO – OFF
- MICROPHONE – OFF
- QUESTIONS – USE CHAT



VIDEO
OFF

MUTE

CHAT





ARCOMA

PREMIUM DIGITAL X-RAY SYSTEMS

AGENDA

- What differs in Aceso and Aceso+ and what do they have in common?
- Overview of system components and how they communicate
- How to navigate our documentation and schematics
- Q&A



WHO AM I?

ERIK SÄLL
ARCOMA ACADEMY MANAGER

- 15 years of experience with medical x-ray equipment
- Hands on technical background (service/installation/troubleshooting)
- Focusing majority of time on training and knowledge transferring activities



WHAT DIFFERS BETWEEN ACESO+ AND ACESO?





ARCOMA


PREMIUM DIGITAL X-RAY SYSTEMS

Aceso+

Arcoma product name = Precision

Most advanced positioning system provided by Arcoma:

- Fully automated system with motorized movements in X, Y, Z, and A, B directions
- Configurable autopositions with various mode functionalities
- Stitching available towards Wallstand and Table
- External service software for calibration




Aceso

Arcoma product name = Intuition

Basic positioning system provided by Arcoma:

- Motorized movements in Z directions only
- Autotracking functionality
- Stitching available towards Wallstand
- Two table types available
- Integrated service software for calibration






CMP200 Generator



CXDI Control Software NE




Overwrap
software




Positioning system with tube screen (Intuition or Precision)

CXDI detectors







Receptor Setup

Receptor

- ☒ R1
- ☐ R2
- ☐ R3
- ☐ R4
- ☐ R5
- ☐ R6

Receptor Properties

General

☒ Receptor Enable

☐ Tomo
 ☐ Serial
 ☐ Fluoro
 ☐ Auto Focus

Memory

☒ Off
 ☐ On
 ☐ Default

Tubes

☒ Tube 1
 ☐ Tube 2

Values

Tomo Back-up Time (ms): 2500

▼

Fluoro Hangover (s): 30

▲

▼

Rad Hangover (s): 0

▲

▼

Last Image Hold (ms): 40

▲

▼

Interface Options: 27

▲

▼

Functional Options: 0

▲

▼

Operations

Apply

Refresh

Exit

Close

Help

Protocol Editor

Protocol

- Protocol
 - A1b Test Wall Stand
 - Wall sim
 - Radiography
 - A1d Test Wall Stand Fault
 - A3 Test Wall Stand
 - A4 Test Wall Stand Coll Fail
 - B1 Test Table Flexible Copy2
 - B1c Test Chest AP Table
 - B1c Test Chest AP Table COF
 - B1d Test Table feltest
 - C Test Universal
 - C1a Test Universal
 - CanonIgnore 1
 - Chest
 - Chest AP C
 - Chest PA
 - Wall sim
 - Radiography
 - Chest PA LAT
 - Chest SID 180cm
 - Chest SID 300cm
 - D1b Test Stitching Wall BetaC
 - D2 Test Stitching Table Copy
 - Demo11 - Stitching Wall 110
 - Demo13 - Shoulder No Wait
 - Demo13 - Stitching Wall, 180
 - Demo14 - Shoulder No Wait
 - Demo3 - Abdomen AP
 - Demo4 - Table Pendulum
 - Demo5 - Test Table stitching
 - Demo6 - Table Flexible Centri
 - Demo7 - Table Flexible Bords

Exposure mode APRID

Radiography kV=40, mA=500, ms=100, Technique=0, Film=0, Focus=0, LeftField=0, CenterField=1, RightField=0, Receptor=3, Density=0, AECFieldsOrientation=0, AutoPosOn=1, AutoPositi

IP Parameter X-ray Parameter

☐ Long exposure

APR-ID: kV=40, mA=500, ms=100, Technique=0, Film=0, Focus=0, LeftField=0, CenterField=1, RightField=0, Receptor=3, Density=0, AECFieldsOrientation=0, AutoPosOn=1, AutoPosit

Fluoro SensorArea MaxPulseWidth SeriesInstanceUID

Binning ADC-ROI

Cine/Ser. Rad. MaxPulseWidth SeriesInstanceUID

Binning ADC-ROI SensorArea

Tomosynthesis Option Tomo Height(mm)

Body Size medium

NAME	Very Small	Small	Medium	Large
Focus	SMALL	SMALL	SMALL	SMALL
Left Field	NO	NO	NO	NO
Center Field	YES	YES	YES	YES
Right Field	NO	NO	NO	NO
Receptor	3	3	3	3
Density	0	0	0	0
AEC Fields Orient.	1-2-3 Portrait	1-2-3 Portrait	1-2-3 Portrait	1-2-3 Portrait
AutoPosition On	YES	YES	YES	YES
Auto Position	12	12	12	12

Which system do you work with? Aceso+ or Aceso?

PLEASE TYPE SYSTEM NAME IN CHAT




Aceso+ (Precision)



Aceso (Intuition)


SYSTEM COMPONENTS *Aceso+*

CAN BUS COMMUNICATION




SYSTEM COMPONENTS *Aceso+*

ETHERNET COMMUNICATION



SYSTEM COMPONENTS *Aceso+*


SERIAL COMMUNICATION



SERVICE TOOL CONNECTIONS


Aceso+

SERIAL COMMUNICATION



SYSTEM COMPONENTS *Aceso+*

PHYSICAL LOCATIONS




OVERVIEW OF SYSTEM COMPONENTS

CLARIFICATIONS

THE TUBE DISPLAY ASSEMBLY CONTAINS
TWO COMMUNICATION MODULES:

- ETHERNET NODE COMMUNICATING WITH CB800

- CAN BUS NODE COMMUNICATING WITH
COLLIMATOR AND Z/MASTER




OVERVIEW OF SYSTEM COMPONENTS

CLARIFICATIONS

THE SYSTEM MASTER NODE (CB800)
COMMUNICATES BOTH BY ETHERNET AND CAN BUS:

- ETHERNET COMMUNICATION WITH TUBE DISPLAY AND CANON PC

- CAN BUS COMMUNICATION WITH Z/MASTER IN OTC




OVERVIEW OF SYSTEM COMPONENTS


COMMUNICATION WITH CANON NE

CPI AND CANON

- CPI Overwrap Software (and also Genware MP – the generator service software) is stored on the Canon PC and communicates with CPI generator
- The CPI Overwrap Software communicates with CANON NE



SERVICE SOFTWARE CONNECTION *Aceso+* *Serial (RS 232)*




CONNECTING SERVICE SOFTWARE: RS 232 COMMUNICATION

When connecting the Arcoma Service Software we use the already wired RS 232 connection on the Z/Master node

CABLE ROUTED FROM
Z/MASTER TO SYSTEM
CABINET FOR EASIER
CONNECTION

CAN bus / Serial RS 232 ports on CBDC-nodes



CONNECTING SERVICE SOFTWARE: RS 232 COMMUNICATION


We must always use RS 232 connection port
for serial communication with any node

NOTE!
Connecting a serial
cable to any CAN
bus port will cause
permanent
damages to node

CAN bus nodes


CONSISTENT ELECTRIC DESIGN:

- DIGITAL I/O = J0, J1, J2, J3
- ANALOGUE I/O = J4
- CAN BUS PORT = J6
- SERIAL PORT = J7
- 24 VDC POWER = J8



CAN BUS SYSTEMS IN GENERAL

CLOSED CIRCUIT




**ALL NODES ON
CAN BUS CIRCUIT
NEEDS TO STAY
CONNECTED**

SOFTWARE UPDATING *DOCUMENTATION*

SOFTWARE PACKAGES AND DOCUMENTATION

- All software packages have a bundle version number

Name	Date modified	Type
 89-109_Canon_CPI (T3)_3.4	11/5/2018 1:11 PM	File folder

**SOFTWARE
DOCUMENTATION:
RLN AND UDI**






SOFTWARE UPDATING *DOCUMENTATION*

RLN AND UDI

All software packages always include two types of critical documents:

- The first one is called RLN
- The second one is called UDI

New folder > 89-109_Canon_CPI (T3)_3.4 > Documentation

Name	Date modified	Type	Size
 SwRLN_0072-C_OTC_2_3.pdf	11/7/2017 3:09 PM	PDF File	507 KB
 SwRLN_0072-C_System_3_4_0.pdf	10/8/2018 1:19 PM	PDF File	535 KB
 SwRLN_H000_1_5_0.pdf	2/18/2014 11:26 AM	PDF File	178 KB
 SwUDI_0072-C_OTC_2_3.pdf	11/8/2017 1:45 PM	PDF File	637 KB
 SwUDI_0072-C_System_3_4.pdf	10/8/2018 1:20 PM	PDF File	1,181 KB

**RLN = RELEASE
NOTES**


**UDI = UPDATE
INSTRUCTION**

SOFTWARE UPDATING *DOCUMENTATION*

RLN - Release notes

Release notes is used to determine which version of each node belonging to a specific software package.

One for OTC-level
(CAN bus positioning system)



2.3 Update Instructions

See document SwUDI_0072-C_OTC_x_x.pdf

3 RELEASE NOTES

3.1 2.3

3.1.1 Software versions

Master/Z:	2.3.A
UI/Collimator:	1.4.A
Guard:	0.10.R
SI:	2.0.A
X:	1.1.B
Y:	1.1.B
Alpha/Beta:	1.0.A
Bucky:	1.3.A
Wallstand:	1.3.A

SSW: 1.5.A

3.1.2 New functions

3.1.3 Enhancements and Bugfixes


- Enhancements in stitching functionality. Overlap and collimator size can be adjusted.

3.1.4 Parameter changes

SOFTWARE UPDATING *DOCUMENTATION*

RLN – Release notes

And another one for system level
(CB800, CPI overwrap)



2.1.2 Software compatibility

0072-C_OTC version 2.3

CANDevice version 1.5

CPI overwrap 3.4.0.0 & 3.6.0.0

2.2 Known Issues

2.3 Update Instructions

See document SwUDI_0072-C_System_x_x.pdf

3 RELEASE NOTES

3.1 3.4

3.1.1 Software versions

SystemMaster: 3.4.A

Display: 3.0.F

3.1.2 New functions

3.1.3 Enhancements and Bugfixes

The handling of preview images is improved.

SOFTWARE UPDATING DOCUMENTATION

UDI – Update Instruction

Same thing with the Update Instructions – there is one for OTC level and another one for system level.

4.1 Update Process

To update the software system:

1. Connect the s
2. Start the servi
3. Under “Edit”
4. Mark all possible setting and save, an example is shown in Figure 1:




Figure 1: Hardware Key Mask

5. Restart the system.

6. Put on power to the board. The diodes on the board shall now follow the sequence:
 - Led 1 starts to blink.
 - Led 1 and 2 starts to blink when USB and file system is detected.
 - Led 1 and 3 starts to blink when the application is programmed.
 - Led 1 and 4 starts to blink when the programming of the board is finished.
7. When the programming is finished, remove the USB-memory stick.
 - The new application starts.




Figure: CB800 ports

SOFTWARE UPDATING

DOCUMENTATION – HOW TO LOAD SW

HOW TO LOAD SOFTWARE

We have created a guidance document explaining the following:

- How to use the RLN and UDI documents
- In which order to read and determine versions
- When and how to save a backup parameter (system config file)
- The difference between upgrading a complete system and replacing a single spare part containing software
- How to load software to both CAN bus and Ethernet nodes

**STEP BY STEP GUIDE
AVAILABLE ON
ARCOMA PARTNER
PORTAL**

SOFTWARE UPDATING

HOW TO LOAD SOFTWARE

Let's have a look at the document together:

[How to load software](#)



**DOWNLOAD
HOW TO
GUIDE**

SCHEMATICS

HOW TO NAVIGATE

There are 3 different types of schematics:

- SBD = System Block Diagram
- UBD = Unit Block Diagram
- WRD = Wiring Diagram


**BLOCK
DIAGRAMS
AND
WIRING
DIAGRAMS**

SCHEMATICS

HOW TO NAVIGATE

All cables and wires are categorized with consistent names:

- 1.x = Tube stand
- 2.x = Table
- 3.x = Wallstand
- 4.x = System cabinet
- 5.x = Image system




SCHEMATICS

HOW TO NAVIGATE EASILY


Different schematic types and consistent labeling
- *Let's have a look at them!*

Installation and Service manual



SCHEMATICS
AND
LABELINGS

Questions?



Feel free to send us questions later:

erik.sall@arcoma.se or service@arcoma.se

THANK YOU

Arcoma AB | Annavägen 1 | 352 46 Växjö | Sweden



PHONE

+46 470 70 69 00



EMAIL

office@arcoma.se



WEBSITE

www.arcoma.se