

# ARCOMA OPEN TABLE



- Outstanding patient coverage with wide detector and table top stroke
- Fully operational at maximum patient weight
- Easy patient loading with low loading position
- High Vertical Speed for fast positioning

# SPECIFICATIONS ARCOMA OPEN TABLE

## GENERAL

- Movement: 6-Way
- Operating method: Motorized vertical and floating table top
- Patient load: 300 kg (660 lbs)
- Motorized vertical travel: 555 - 930 mm (22" - 36")
- Longitudinal travel:  $\pm 600$  mm ( $\pm 23$ " )
- Lateral travel:  $\pm 150$  mm ( $\pm 6$ " )

## TABLE TOP HEIGHT:

- Lowest table top position (from floor to table top surface): 55 cm
- Vertical stroke: 38 cm

## CONTROLS

- Table hand control
- Foot switch (releases long. and lat. separately)
- Foot switch (Floor placed 1260mm long table top release switch)

## OPTION/ACCESSORIES:

- Vertical collision protection
- Compression belt
- Patient hand grip
- Mattress