



WORLD CLASS DIGITAL X-RAY SYSTEMS

LAUNCH PRESENTATION

ARCOMA PRECISION i5

WE ARE HAPPY TO PRESENT THE NEW:

ARCOMA PRECISION i5



Powered by
Canon
DR

THE NEW PRODUCT

Arcoma are happy to launch its new product Arcoma Precision i5. The new product consists of a brand new wall stand with motorized tilt for better ergonomics and optimal workflow together with a new OTC with brand new modern Scandinavian design and smart features for easy and quick use.

The new Precision i5 is also built on a silent concept for optimal patient experience and optimal working environment.





Motorized Wall stand

+

New OTC

+

Silent Concept

+

Powered by
Canon
DR

MOTORIZED WALL STAND

Optimized workflow

Motorized movements enables efficient examinations

Ergonomic design

Motorized tiltable wall stand allows minimal manual handling.

Increased focus on patient

Easy adjustment and motorized tilt enables increased focus on the patient.




NEW MOTORIZED WALL STAND

BENEFITS	FEATURES
<ul style="list-style-type: none">• Higher patient throughput thanks to faster positioning• Highest and lowest detector position on the market• Lighter lateral armrest with adjustable angle• Faster installation	<ul style="list-style-type: none">• Motorized tilt of the detector• Motorized and manual vertical movements• Auto tracking vertical and tilting movements• Pendulum• Hand control for operating wall stand, collimator and auto positioning• Removable lateral arm rest with parking position and adjustable angle



NEW MOTORIZED WALL STAND

EASY ACCESSIBLE INTERFACE



NEW OTC DESIGN

Excellent design

Modern design that blends in nicely into any hospital environment.

Easy disinfection

Easy cleaning and disinfection thanks to glass overlay display with no edges.

Smart and easy to use

Good viewing angle of display and few buttons makes the OTC user-friendly and easy to understand.



NEW OTC DESIGN

BENEFITS	FEATURES
<ul style="list-style-type: none">• Simplified work - easier to use• Blends nicely into any hospital room• All data displayed in one display• Modern user interface• Better ergonomics with new buttons• Easy access buttons• Tactile recognition buttons for safe and fast operation	<ul style="list-style-type: none">• New design• Easy cleaning• Integrated light indication• All directions movement• Adapted for left- and righthanded




NEW OTC DESIGN

FLEXIBLE USER INTERFACES

ALL DATA IN ONE INTERFACE

Release in all directions
- easy movements



Up / down easy accessible
for both left- and righthanded

Tactile
recognition

Light indication

Release
- all directions

Up / down

EASY ACCESSIBLE INTERFACE



FULL INTEGRATION

- Patient information
- Exposure parameters
- Collimator parameters

Release
- all directions
Release in
one direction



SILENT CONCEPT

Improved patient experience and user work environment

- Minimized sounds for better work environment and better patient experience
- New power supplies without fans in both table, OTC and wall stand



NEW SILENT CONCEPT

BENEFITS	FEATURES
<ul style="list-style-type: none">• Improved patient experience• Improved patient environment• Improved working environment	<ul style="list-style-type: none">• New power supplies without fans in OTC, table and wall stand



THANK YOU

Arcoma AB | Annavägen 1 | 352 46 Växjö | Sweden



PHONE

+46 470 70 69 00



EMAIL

office@arcoma.se



WEBSITE

www.arcoma.se